

Labycare-Ultrasound Therapy Unit Digital Dual Head (SKSS-PP341)



Ultrasound Therapy Unit Digital Dual Head

Description :

The microprocessor based, 1MHZ dual frequency unit with LCD Display & a Soft touch membrane keyboard enables the physician to set the parameters easily.

Technical :

Output Frequency : Two, 1MHz and 3MHz

Output Watts : Continuous Mode : 15 watts, Pulsed Mode : 21 watts

Intensity : 0 to 3.5 watts / cm² maximum. US output can be adjusted by intensity switch.

Treatment Heads : Two. (One Small & One Large of 5cm sq.)

Pulse Mode Ratios : 1:1 , 1:2 , 1:4 , 1:7 , 1:10 , 1:16

Effective Area of Radiation : 5 square centimeter (Large Head).

Treatment Timer : 00 to 60 minutes with Buzzer at the end of Treatment Time.

Display : The LCD displays Parameter chosen, Treatment Time and Output Power.

Programs : 45 pre-fixed programmed

Power : 230V / AC / 50Hz. 0.6Amps.

We are Exporters and Suppliers of Ultrasound Therapy Unit Digital Dual Head at wholesale competitive prices from India.

Specifications :

The microprocessor based, 1MHZ dual frequency unit with LCD Display & a Soft touch membrane keyboard enables the physician to set the parameters easily.

Technical :

Output Frequency : Two, 1MHz and 3MHz

Output Watts : Continuous Mode : 15 watts, Pulsed Mode : 21 watts

Intensity : 0 to 3.5 watts / cm² maximum. US output can be adjusted by intensity switch.

Treatment Heads : Two. (One Small & One Large of 5cm sq.)

Pulse Mode Ratios : 1:1 , 1:2 , 1:4 , 1:7 , 1:10 , 1:16

Effective Area of Radiation : 5 square centimeter (Large Head).

Treatment Timer : 00 to 60 minutes with Buzzer at the end of Treatment Time.

Display : The LCD displays Parameter chosen, Treatment Time and Output Power.

Programs : 45 pre-fixed programmed

Power : 230V / AC / 50Hz. 0.6Amps.

We are Exporters and Suppliers of Ultrasound Therapy Unit Digital Dual Head at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com