

Labycare-Interferential Therapy Unit (SKSS-PP343)



Interferential Therapy Unit MANUFACTURER BY LABYCARE

Description :

Interferential current stimulation is effective upon the Nerves, Muscles, Arteries and Veins.

The current is able to penetrate tissues with a soft stimulative effect. thus assuring an exclusive range of application.

Technical :

MEDIUM FREQUENCY : 2 KHz and 4KHz & Russian Curren.

BASE AMF : 0 to 100Hz

SWEEP AMF : 0 to 15Hz

SPECTRUM (From) : 3 Position : Rectangular 1/1 Sec, Trapezoidal 1/5 Sec, Triangular 6/6 Sec.

MUSCLE STIMULATOR : One Channel

THERAPY MODE :

(a) Four Pole, 45 deg. Vector IFT.

(b) Four Pole, 90 deg. Vector IFT

(c) Four Pole IFT

(d) Bi Pole IFT (2 Electrodes)

INTENSITY : 0-60mA Continuously Adj.

INTENSITY BALANCE : Available Equalize current between Two Channels.

TIMER : 00 to 60 minutes, Digital

POWER SUPPLY : 230V AC, 50Hz

BUZZER : audio buzzer at the end of Treatment Time.

Accessories :

Gel one bottle,

Velcro Straps to hold / tie Pads-2 , User Manual.

We are Exporters and Suppliers of Interferential Therapy Unit at wholesale competitive prices from India.

Specifications :

Interderential current stimulation is effective upon the Nerves, Muscles, Arteries and Veins.

The current is able to penetrate tissues with a soft stimulative effect. thus assuring an exclusive range of application.

Technical :

MEDIUM FREQUENCY : 2 KHz and 4KHz & Russian Curren.

BASE AMF : 0 to 100Hz

SWEEP AMF : 0 to 15Hz

SPECTRUM (From) : 3 Position : Rectangular 1/1 Sec, Trapezoidal 1/5 Sec, Triangular 6/6 Sec.

MUSCLE STIMULATOR : One Channel

THERAPY MODE :

(a) Four Pole, 45 deg. Vector IFT.

(b) Four Pole, 90 deg. Vector IFT

(c) Four Pole IFT

(d) Bi Pole IFT (2 Electrodes)

INTENSITY : 0-60mA Continuously Adj.

INTENSITY BALANCE : Available Equalize current between Two Channels.

TIMER : 00 to 60 minutes, Digital

POWER SUPPLY : 230V AC, 50Hz

BUZZER : audio buzzer at the end of Treatment Time.

Accessories :

Gel one bottle,

Velcro Straps to hold / tie Pads-2 , User Manual.

We are Exporters and Suppliers of Interferential Therapy Unit at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com