

Labycare-Physio Suspension Frame without Couch (SKSS-PP349)



Physio Suspension Frame without Couch Manufacturer by LABYCARE

Description :

Tubular Steel Construction, with welded steel rod-mesh at the top for complete suspension with slings & ropes.

Mesh facilitates suspension at any desired point in relation to the patient and as required in each individual case. Comes complete with suspension Gear.

Frame supplied without couch.

Technical :

Construction : Tubular Steel Frame with 6mm thick rod-mesh at the top and on upper sides. Legs fitted with Rubber / PVC Shoes.

Size : Overall sizes 210cm long X 110Ccm wide at the base X 50cm wide at the top & 190cm high.

Suspension Gear : Comprising 7 canvas slings(4 limbs, 2 Trunks, & 1 Head), 4 Single Hook Ropes, 2 Single Pulley Ropes,& 2 Double Pulley Ropes. with Necessary Hooks & Cleats, 8 Graduated Springs of Different Marked Tension. 2 Ankles Straps and 4 'S' Hooks.

Finish : Off white / Cream, Powder Coated Finish.

We are Exporters and Suppliers of Physio Suspension Frame without Couch at wholesale competitive prices from India.

Specifications :

Tubular Steel Construction, with welded steel rod-mesh at the top for complete suspension with slings

& ropes.

Mesh facilitates suspension at any desired point in relation to the patient and as required in each individual case. Comes complete with suspension Gear.

Frame supplied without couch.

Technical :

Construction : Tubular Steel Frame with 6mm thick rod-mesh at the top and on upper sides. Legs fitted with Rubber / PVC Shoes.

Size : Overall sizes 210cm long X 110Ccm wide at the base X 50cm wide at the top & 190cm high.

Suspension Gear : Comprising 7 canvas slings(4 limbs, 2 Trunks, & 1 Head), 4 Single Hook Ropes, 2 Single Pulley Ropes,& 2 Double Pulley Ropes. with Necessary Hooks & Cleats, 8 Graduated Springs of Different Marked Tension. 2 Ankles Straps and 4 'S' Hooks.

Finish : Off white / Cream, Powder Coated Finish.

We are Exporters and Suppliers of Physio Suspension Frame without Couch at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com