

Labycare-Suspension Spring with Marked (SKSS-PP351)



Suspension Spring with Marked Manufacturer by LABYCARE

Description :

Springs : Suspension Springs are made out of best steel wire, with chrome plated finish, rope is fitted in the spring to avoid over stretching.

Tension of the springs are about : A,5kg-7kg, B,8kg-10kg, C,12kg-15kg, D,18kg-20kg, E,22kg-25kg.

Set consists of ONE spring each (A+B+C+D+E) with "S" hooks.

Suitable for extension, Flexion and Suspension Activities.

Full Set of 5 Springs.

We are Exporters and Suppliers of Suspension Spring with Marked at wholesale competitive prices from India.

Specifications :

Springs : Suspension Springs are made out of best steel wire, with chrome plated finish, rope is fitted in the spring to avoid over stretching.

Tension of the springs are about : A,5kg-7kg, B,8kg-10kg, C,12kg-15kg, D,18kg-20kg, E,22kg-25kg.

Set consists of ONE spring each (A+B+C+D+E) with "S" hooks.

Suitable for extension, Flexion and Suspension Activities.

Full Set of 5 Springs.

We are Exporters and Suppliers of Suspension Spring with Marked at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com