

Labycare-Exercising Unit Multi Therapy (SKSS-PP353)



Exercising Unit Multi Therapy Manufacturer & Supplier

Description :

Bright Polished Hardwood constructed Multi-Exercise unit with pulley system is suitable for shoulder, Elbow, Wrist, Chest and Lower limb exercises.

Unit Consists of the following attachments :

- (a) Wall Bars 240cm high x 75cm wide with knee support pad.
- (b) Detachable wrist roller with side handle for supination and pronation exercises.
- (c) Removable Rowing Machine cum sliding seat.
- (d) Unit is fitted with 12 pulleys and Two sets of weight plates (10kg in each set) that allows effective weight resistive exercises at different levels.

We are Exporters and Suppliers of Exercising Unit Multi Therapy at wholesale competitive prices from India.

Specifications :

Bright Polished Hardwood constructed Multi-Exercise unit with pulley system is suitable for shoulder, Elbow, Wrist, Chest and Lower limb exercises.

Unit Consists of the following attachments :

- (a) Wall Bars 240cm high x 75cm wide with knee support pad.
- (b) Detachable wrist roller with side handle for supination and pronation exercises.
- (c) Removable Rowing Machine cum sliding seat.

(d) Unit is fitted with 12 pulleys and Two sets of weigh plates (10kg in each set) that allows effective weight resistive exercises at different levels.

We are Exporters and Suppliers of Exercising Unit Multi Therapy at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com