

## Labycare-Arm Exercise Unit (SKSS-PP356)



Arm Exercise Unit Manufacturer & Supplier

### **Description :**

Arm Exerciser that allows effective weight resistive exercise at different levels for shoulder, Elbow and Chest.

### **Technical :**

**Construction :** Chrome Plated, Tubular Steel construction mounted on wooden platform and fitted with painted wall-bracket at the top for wall mounting.

**Weight Plates :** One Set (10 X 1kg plate in set)

**Handles :** Two, (One for single grip & One for both Hands).

**Pulley Block :** Swivel pulley block with Height Adjustment at different levels.

**Rope Adjuster :** Rope length can be adjusted as per Height level of Pulley Block.

**Ankle Strap :** One.

**Wall Mounting :** Unit Comes with wall fixing hardware.

We are Exporters and Suppliers of Arm Exercise Unit (Pulley Apparatus) at wholesale competitive prices from India.

### **Specifications :**

Arm Exerciser that allows effective weight resistive exercise at different levels for shoulder, Elbow and Chest.

### **Technical :**

**Construction** : Chrome Plated, Tubular Steel construction mounted on wooden platform and fitted with painted wall-bracket at the top for wall mounting.

**Weight Plates** : One Set (10 X 1kg plate in set)

**Handles** : Two, (One for single grip & One for both Hands).

**Pulley Block** : Swivel pulley block with Height Adjustment at different levels.

**Rope Adjuster** : Rope length can be adjusted as per Height level of Pulley Block.

**Ankle Strap** : One.

**Wall Mounting** : Unit Comes with wall fixing hardware.

We are Exporters and Suppliers of Arm Exercise Unit (Pulley Apparatus) at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)