

Labycare-Wall Bar Wooden (SKSS-PP357)



Wall Bar Wooden Manufacturer by LABYCARE

Description :

Wall bar single section, hardwood constructed with sturdy frame finished in natural wood polish.

Oval shaped rungs are squarely fitted into vertical bars. Patient can perform shoulder range of motion exercises from sitting or standing position.

Unit is provided with two wall boards for fixing.

Technical :

Size : 90cm wide X 244cm high

Rungs : 13, oval shaped rungs fitted squarely into vertical bars.

Adjustable Bars : One Bar (10th) is adjustable in height in Three steps.

Finish : Natural wood polish finish.

Chinning Bar : One bar is provided at top outside the frame for chinning exercise.

Standard Accessories :

Knee Support : One, foam padded knee support is provided as an aid for wheelchair patients for pull up exercises.

We are Exporters and Suppliers of Wall Bar Wooden at wholesale competitive prices from India.

Specifications :

Wall bar single section, hardwood constructed with sturdy frame finished in natural wood polish.

Oval shaped rungs are squarely fitted into vertical bars. Patient can perform shoulder range of motion exercises from sitting or standing position.

Unit is provided with two wall boards for fixing.

Technical :

Size : 90cm wide X 244cm high

Rungs : 13, oval shaped rungs fitted squarely into vertical bars.

Adjustable Bars : One Bar (10th) is adjustable in height in Three steps.

Finish : Natural wood polish finish.

Chinning Bar : One bar is provided at top outside the frame for chinning exercise.

Standard Accessories :

Knee Support : One, foam padded knee support is provided as an aid for wheelchair patients for pull up exercises.

We are Exporters and Suppliers of Wall Bar Wooden at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com