

Labycare-Wall Bar cum Multi Exercise Unit (SKSS-PP358)



Wall Bar cum Multi Exercise Unit Manufacturer by LABYCARE

Description :

Multi exerciser cum wall bar unit is in single section & is made of tubular pipe with powder coated finish. Patient can perform shoulder range of motion exercises from sitting or standing position.

Unit comes with attachment for Pull-up exercises, shoulder Pulley attachment & shoulder abduction ladder.

Technical :

Construction : Rectangular and Round Pipe constructed unit with hooks for shoulder pulleys & plastic steps for shoulder abduction exercises.

Rungs : 13, multi-colour, round metal rungs are fitted into vertical rectangular pipes.

Chinning Bar : One bar is provided at top outside the frame for chinning exercises.

Knee Support : One, Foam padded knee support is provided as an aid for wheelchair patients for pull-up exercises.

Shoulder Pulleys : A Set of pulley with handles are provided for over-head pulley exercises.

Shoulder Abduction : 30, Plastic Finger holds in different colours are provided for shoulder abduction exercises.

Size : 90cm wide X 240cm high.

Unit is provided with two wall plates with wall fixing hardware.

Finish : Powder Coated Finish.

We are Exporters and Suppliers of Wall Bar cum Multi Exercise Unit at wholesale competitive prices from India.

Specifications :

Multi exerciser cum wall bar unit is in single section & is made of tubular pipe with powder coated finish. Patient can perform shoulder range of motion exercises from sitting or standing position.

Unit comes with attachment for Pull-up exercises, shoulder Pulley attachment & shoulder abduction ladder.

Technical :

Construction : Rectangular and Round Pipe constructed unit with hooks for shoulder pulleys & plastic steps for shoulder abduction exercises.

Rungs : 13, multi-colour, round metal rungs are fitted into vertical rectangular pipes.

Chinning Bar : One bar is provided at top outside the frame for chinning exercises.

Knee Support : One, Foam padded knee support is provided as an aid for wheelchair patients for pull-up exercises.

Shoulder Pulleys : A Set of pulley with handles are provided for over-head pulley exercises.

Shoulder Abduction : 30, Plastic Finger holds in different colours are provided for shoulder abduction exercises.

Size : 90cm wide X 240cm high.

Unit is provided with two wall plates with wall fixing hardware.

Finish : Powder Coated Finish.

We are Exporters and Suppliers of Wall Bar cum Multi Exercise Unit at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com