

Labycare-Hospital Examination Table (SKSS864A)



Hospital Examination Table manufacturer by **Labycare** in India.

Description : **Labycare Hospital Examination Table** is used to support patients during medical examinations. During these exams, doctors in offices, clinics and hospitals use an adjusting mechanism to manipulate and position the table to allow patient support, closer examination of a portion or the entire patient, and the ability to move the patient on and off the table safely.

Specifications :

- Back section adjustable to 75 degree
- Powder coated steel frame, firm and anti-rust
- With a shelf to assemble non woven, easy to clean and more anti-infection
- High quantity leather, it is waterproof
- Size: 190 x 80 x 45 cm
- Other size and material offered

labycare

www.labycare.com | Email: labycare@gmail.com