

Labycare-Shoulder Double Pulley Set (SKSS-PP362)



Shoulder Double Pulley Set Manufacturer by LABYCARE

Description :

Tubular Steel constructed wall mounting unit provides reciprocal, assistive & resistive exercises of the elbow & shoulder.

Technical :

Construction : Tubular Steel bracket of 20mm dia. tube is fitted with Two close hooks to fix pulleys.

Pulleys : Two, Nylon pulleys covered in metal sheet casing with Hooks to fix..

Grip Handles : Two, spring steel wire handles with plastic grips.

Rope / Cord : Durable Interwoven Nylon Cord of Suitable length.

Mounting : Wall mountable unit, set include mounting hardware.

Finish : Powder Coated Finish.

We are Exporters and Suppliers of Shoulder Double Pulley Set at wholesale competitive prices from India.

Specifications :

Tubular Steel constructed wall mounting unit provides reciprocal, assistive & resistive exercises of the elbow & shoulder.

Technical :

Construction : Tubular Steel bracket of 20mm dia. tube is fitted with Two close hooks to fix pulleys.

Pulleys : Two, Nylon pulleys covered in metal sheet casing with Hooks to fix..

Grip Handles : Two, spring steel wire handles with plastic grips.

Rope / Cord : Durable Interwoven Nylon Cord of Suitable length.

Mounting : Wall mountable unit, set include mounting hardware.

Finish : Powder Coated Finish.

We are Exporters and Suppliers of Shoulder Double Pulley Set at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com