

Labycare-shoulder abduction ladder (SKSS-PP363)



shoulder abduction ladder labycare

Description :

Shoulder and Finger Ladder made of polished Teak-wood having wooden steps for shoulder abduction exercises.

Unit Size : 8cm wide X 6cm deep X 100cm high.

Steps : Thirty wooden steps & each step is numbered to keep record of progress.

Mounting : Wall mountable unit, set include mounting hardware.

Finish : Natural wood polish finish.

We are Exporters and Suppliers of Shoulder Abduction Ladder at wholesale competitive prices from India.

Specifications :

Shoulder and Finger Ladder made of polished Teak-wood having wooden steps for shoulder abduction exercises.

Unit Size : 8cm wide X 6cm deep X 100cm high.

Steps : Thirty wooden steps & each step is numbered to keep record of progress.

Mounting : Wall mountable unit, set include mounting hardware.

Finish : Natural wood polish finish.

We are Exporters and Suppliers of Shoulder Abduction Ladder at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com