

Labycare-Parallel Walking Bars Adult (SKSS-PP366)



Parallel Walking Bars Adult Manufacturer By LABYCARE

Description :

Length 3.60 meters, with Hand rails made of PVC sheathed iron tubing.

Height & Width adjusting arms made of chrome plated steel tubing of 38 mm dia.

Outer uprights of thick steel tubing (50mm dia).

Mounted on hard wood platform with inclination on both ends and covered with non-slip matting.

Height is adjustable from 76mm to 10cm and width from 40cm to 65 cm.

Adjustment can easily be done by means of setting pins at the required holes and can further be stabilized by tightening the screws.

Comes with detachable polished hard wood abduction board. Steel uprights are oven baked finish.

We are Exporters and Suppliers of Parallel Walking Bars Adult at wholesale competitive prices from India.

Specifications :

Length 3.60 meters, with Hand rails made of PVC sheathed iron tubing.

Height & Width adjusting arms made of chrome plated steel tubing of 38 mm dia.

Outer uprights of thick steel tubing (50mm dia).

Mounted on hard wood platform with inclination on both ends and covered with non-slip matting.

Height is adjustable from 76mm to 10cm and width from 40cm to 65 cm.

Adjustment can easily be done by means of setting pins at the required holes and can further be stabilized by tightening the screws.

Comes with detachable polished hard wood abduction board. Steel uprights are oven baked finish.

We are Exporters and Suppliers of Parallel Walking Bars Adult at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com