

Labycare-Parallel Walking Bars Child (SKSS-PP367)



Parallel Walking Bars Child Manufacturer by LABYCARE

Description :

Hand Rails made of PVC sheathed iron tubing.

Mounted on a 3 Meter long wooden platform with inclination in both ends and covered with non-slip matting. Height and width adjustments can easily be done by means of setting pins & screws.

Comes with detachable polished wooden abduction board.

Height is adjustable from 53cm to 82cm and width from 35cm to 54cm.

Adjustment can easily be done by means of setting pins the required holes and can further be stabilized by tightening the screws.

We are Exporters and Suppliers of Parallel Walking Bars Child at wholesale competitive prices from India.

Specifications :

Hand Rails made of PVC sheathed iron tubing.

Mounted on a 3 Meter long wooden platform with inclination in both ends and covered with non-slip matting. Height and width adjustments can easily be done by means of setting pins & screws.

Comes with detachable polished wooden abduction board.

Height is adjustable from 53cm to 82cm and width from 35cm to 54cm.

Adjustment can easily be done by means of setting pins the required holes and can further be stabilized by tightening the screws.

We are Exporters and Suppliers of Parallel Walking Bars Child at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com