

Labycare-Parallel Bar Adult without Platform (SKSS-PP368)



Parallel Bar Adult without Platform Manufacturer by LABYCARE

Description :

Length 3 meters, with bars Hand rails made of PVC sheathed iron tubing.

Height adjusting arms made of chrome plated steel tubing.

Height is adjustable from 76cm to 100cm.

Uprights are mounted on two non-slip aluminium bases with iron angle reinforcement.

Adjustment can easily be done by means of setting pins at the required holes and can further be stabilized by tightening the screws.

Width of the bars are 660mm . Upright finished oven baked in blue.

We are Exporters and Suppliers of Parallel Bar Adult without Platform at wholesale competitive prices from India.

Specifications :

Length 3 meters, with bars Hand rails made of PVC sheathed iron tubing.

Height adjusting arms made of chrome plated steel tubing.

Height is adjustable from 76cm to 100cm.

Uprights are mounted on two non-slip aluminium bases with iron angle reinforcement.

Adjustment can easily be done by means of setting pins at the required holes and can further be stabilized by tightening the screws.

Width of the bars are 660mm . Upright finished oven baked in blue.

We are Exporters and Suppliers of Parallel Bar Adult without Platform at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com