

Labycare-Parallel Walking Bar Deluxe (SKSS-PP369)



Parallel Walking Bar Deluxe Manufacturer by LABYCARE

Description :

Specifically designed deluxe Parallel Bar provides independent Height Adjustment with a touch of a Button. The quick, smooth, electric powered movement to different heights helps to achieve maximum comfort for each patient & a real time saver to the therapists.

Glider provides assistance to patients with hemiplegia, xtremity or trunk weakness, balance problems or those who cannot bear weight on their wrist or hands, width manually adjustable.

Technical :

Hand Rails : Stainless Steel 3.60 meter Long & of 38mm dia.

Height Adjustment : From 75cm to 100cm. Electrically 220v AC.

Width Adjustment : Manually from 40cm to 60cm.

Glider Parallel Bar : Padded Forearm support with grip handles, Positioned and moves a top bars.

Abduction Board : Wooden Removable with natural finish.

Platform (Base) : Hardwood construction with Tapered ends & fitted with non-slip mattress.

Distance : A Scale in provided on both side to measure & second the distance traveled.

We are Exporters and Suppliers of Parallel Walking Bar Deluxe at wholesale competitive prices from India.

Specifications :

Specifically designed deluxe Parallel Bar provides independent Height Adjustment with a touch of a

Button. The quick, smooth, electric powered movement to different heights helps to achieve maximum comfort for each patient & a real time saver to the therapists.

Glider provides assistance to patients with hemiplegia, xtremity or trunk weakness, balance problems or those who cannot bear weight on their wrist or hands, width manually adjustable.

Technical :

Hand Rails : Stainless Steel 3.60 meter Long & of 38mm dia.

Height Adjustment : From 75cm to 100cm. Electrically 220v AC.

Width Adjustment : Manually from 40cm to 60cm.

Glider Parallel Bar : Padded Forearm support with grip handles, Positioned and moves a top bars.

Abduction Board : Wooden Removable with natural finish.

Platform (Base) : Hardwood construction with Tapered ends & fitted with non-slip mattress.

Distance : A Scale in provided on both side to measure & second the distance traveled.

We are Exporters and Suppliers of Parallel Walking Bar Deluxe at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com