

Labycare-HI-LO ICU BED (SKSS801A2)



HI-LO ICU BED manufacturer by Labycare in India. As a quality oriented organization, we make sure to manufacturer **Fully Electrical Hospital Beds** with central braking and handset remote control system. These **Automatic Hospital Beds** are adjustable according to the comfort of patients.

Description :

Labycare HI-LO ICU BED

Backed by the team of expert professionals we are engaged in manufacturing, supplying and exporting ICU Bed. Offered beds are known for their robust construction, comfort level, adjustable height and durability. These beds are made in strict compliance with the international standards. Once our ICU Beds are manufactured, these are tested on various parameters under the supervision of our expert team. Customers can avail these beds from us at market leading prices.

Specifications :

- Four pieces individual side-rails, ABS moulded
- Heavy ABS Head and Foot Panel with locking
- 4 section top made up of perforated CRCA M.S. sheet
- Framework made up of M.S. tubes
- Bed mounted on 5" ABS Castors, 2 with brakes
- Holder for Urine Bag
- Sitting position adjustable
- Finish: Pre-treated & Epoxy Powder Coated
- Overall approx. size: 2090L x 900W x 580-750H mm
- Positions obtained by Smooth Crank System with foldable handle:

1. Hi-Low

2. Backrest elevation

3. Kneerest elevation

4. Trendelenburg/ Reverse Trendelenburg

Optional

- C4 Section Mattress
- S.S. Telescopic I.V. Rod with 4 location
- Stainless Steel/ Mild Steel Railing
- Stainless Steel Head & Foot Bows

labycare

www.labycare.com | Email: labycare@gmail.com