

Labycare-rehab trainer (SKSS-PP373)



rehab trainer

Description :

This Unique Upper and Lower extremity cycle accommodate users with varying abilities.

Based on German technology it allows the user to passively work their lower body while actively working their upper body or vice - versa.

Technical :

Construction : Heavy Steel frame for exceptional stability & fitted with heavy duty Ball-bearing for easy operation .

Angle Adjustment : Working Angle of the cycle can be adjusted from 45 to 85 degree for comfort of patient.

Seat : Padded seat is adjustable in Height & Distance to accommodate small or tall adults.

Back Support : Padded Back support is adjustable in height & Distance for patient's comfort.

Resistance : A control Knob adjusts resistance with very low start up load.

Display : A fitness computer display and monitors Time, Distance, Speed and Calories.

Foot Hold : Foot Hold Paddles with straps are provided to hold the foot.

Accessories :

Unit Comes with one set of hand holding mitts and foot hold plates.

We are Exporters and Suppliers of Rehab Trainer at wholesale competitive prices from India.

Specifications :

This Unique Upper and Lower extremity cycle accommodate users with varying abilities.

Based on German technology it allows the user to passively work their lower body while actively working their upper body or vice - versa.

Technical :

Construction : Heavy Steel frame for exceptional stability & fitted with heavy duty Ball-bearing for easy operation .

Angle Adjustment : Working Angle of the cycle can be adjusted from 45 to 85 degree for comfort of patient.

Seat : Padded seat is adjustable in Height & Distance to accommodate small or tall adults.

Back Support : Padded Back support is adjustable in height & Distance for patient's comfort.

Resistance : A control Knob adjusts resistance with very low start up load.

Display : A fitness computer display and monitors Time, Distance, Speed and Calories.

Foot Hold : Foot Hold Paddles with straps are provided to hold the foot.

Accessories :

Unit Comes with one set of hand holding mitts and foot hold plates.

We are Exporters and Suppliers of Rehab Trainer at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com