

Labycare-Training Stair Single Side (SKSS-PP378)



Training Stair Single Side Manufacturer By LABYCARE

Description :

Sturdily built in two sections for installation in a straight single side.

Steps arrangement :

Four 6" (15cm) high steps lead to a Platform size 30" x 30" 24" high.

Steps : 30" wide and 11" deep.

Hand Rails : All around hardwood handrails at two different heights to accomodate Adults and Children.

Construction : Steps and platform made of commercial board, hand rails ith supporting bars made out of hardwood. Polished natural wood finish.

We are Exporters and Suppliers of Training Stair Single Side at wholesale competitive prices from India.

Specifications :

Sturdily built in two sections for installation in a straight single side.

Steps arrangement :

Four 6" (15cm) high steps lead to a Platform size 30" x 30" 24" high.

Steps : 30" wide and 11" deep.

Hand Rails : All around hardwood handrails at two different heights to accommodate Adults and Children.

Construction : Steps and platform made of commercial board, hand rails with supporting bars made out of hardwood. Polished natural wood finish.

We are Exporters and Suppliers of Training Stair Single Side at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com