

Labycare-Exercise Staircase Straight Type 75cm (SKSS-PP380)



Exercise Staircase Straight Type 75cm Manufacturer by LABYCARE

Description :

Sturdily Built in two section for Installation in a Straight Line.

Steps arrangement.

Four 15cm high steps on one side leading to a platform of 75cm x 75cm.

(First step moves into the second to make it bus step of 30cm height).

Eight 7.5cm steps on the other side.

Steps are 75cm wide & 22.5cm deep.

Two sets of Hardwood handrails at different height to accommodate Adult and Children.

Step and Platform made of commercial board of 12mm thickness is covered with non-slip matting.

Sides and supporting frame of steps is made of commercial board of 18mm thickness.

Hand Rails & Supporting bars are made of hardwood, Polished Natural Wood Finish.

We are Exporters and Suppliers of Exercise Staircase Straight Type 75cm Wide at wholesale competitive prices from India.

Specifications :

Sturdily Built in two section for Installation in a Straight Line.

Steps arrangement.

Four 15cm high steps on one side leading to a platform of 75cm x 75cm.

(First step moves into the second to make it bus step of 30cm height).

Eight 7.5cm steps on the other side.

Steps are 75cm wide & 22.5cm deep.

Two sets of Hardwood handrails at different height to accommodate Adult and Children.

Step and Platform made of commercial board of 12mm thickness is covered with non-slip matting.

Sides and supporting frame of steps is made of commercial board of 18mm thickness.

Hand Rails & Supporting bars are made of hardwood, Polished Natural Wood Finish.

We are Exporters and Suppliers of Exercise Staircase Straight Type 75cm Wide at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com