

Labycare-Exercise Staircase Straight Type with Ramp 60cm (SKSS-PP382)



Exercise Staircase Straight Type with Ramp 60cm Manufacturer By Labycare

Description :

Sturdily Built in two section for Installation in a Straight Line.

Steps arrangement.

One Section has five steps i.e. 7.5cm, 10cm, 12.5cm, 15cm, & 17.5cm in height and 30cm deep with platform of 60cm x 60cm.

The Other section has a gradual slope 180cm long.

Two sets of hardwood hand rails at different height to accommodate adult and children.

Steps and platform made of commercial board of 12mm thickness is covered with non-slip matting.

Side and supporting frame of steps made of commercial board of 18mm thickness.

Hand rails and Supporting bars are made of hardwood. Polished natural wood finish.

We are Exporters and Suppliers of Exercise Staircase Straight Type with Ramp 60cm Wide at wholesale competitive prices from India.

Specifications :

Sturdily Built in two section for Installation in a Straight Line.

Steps arrangement.

One Section has five steps i.e. 7.5cm, 10cm, 12.5cm, 15cm, & 17.5cm in height and 30cm deep with platform of 60cm x 60cm.

The Other section has a gradual slope 180cm long.

Two sets of hardwood hand rails at different height to accommodate adult and children.

Steps and platform made of commercial board of 12mm thickness is covered with non-slip matting.

Side and supporting frame of steps made of commercial board of 18mm thickness.

Hand rails and Supporting bars are made of hardwood. Polished natural wood finish.

We are Exporters and Suppliers of Exercise Staircase Straight Type with Ramp 60cm Wide at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com