

Labycare-Examination Couch Wooden (SKSS-PP386)



Examination Couch Wooden Manufacturer by LABYCARE

Description :

Standard Treatment cum Examination table with vinyl upholstered well padded top, Size : 72" L x 24" W x 31" H, is provided with adjustable backrest of 18" length which can be adjusted in height for desired position.

One writing board, one drawer with small cub board is provided underneath the top.

Top is mounted on strudily built, natural finished hardwood solid legs.

We are Exporters and Suppliers of Examination Couch Wooden at wholesale competitive prices from India.

Specifications :

Standard Treatment cum Examination table with vinyl upholstered well padded top, Size : 72" L x 24" W x 31" H, is provided with adjustable backrest of 18" length which can be adjusted in height for desired position.

One writing board, one drawer with small cub board is provided underneath the top.

Top is mounted on strudily built, natural finished hardwood solid legs.

We are Exporters and Suppliers of Examination Couch Wooden at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com