

Labycare-Treatment Therapy Couch Wooden (SKSS-PP387)



Treatment Therapy Couch Wooden Manufacturer by LABYCARE

Description :

Standard Therapy Couch with vinyl upholstered single section top.

Technical :

Construction : Top is mounted on strudily built Teak-wood solid legs with H-brace to provide stable treatment surface.

Top : Top of 19mm thick ply with upholstered 50mm thick well padded top & is covered with good quality rexine.

Size : Top 200cm long X 70cm wide X 80cm high.

Finish : Natural Wood Polish Finish.

We are Exporters and Suppliers of Treatment Therapy Couch Wooden at wholesale competitive prices from India.

Specifications :

Standard Therapy Couch with vinyl upholstered single section top.

Technical :

Construction : Top is mounted on strudily built Teak-wood solid legs with H-brace to provide stable treatment surface.

Top : Top of 19mm thick ply with upholstered 50mm thick well padded top & is covered with good quality rexine.

Size : Top 200cm long X 70cm wide X 80cm high.

Finish : Natural Wood Polish Finish.

We are Exporters and Suppliers of Treatment Therapy Couch Wooden at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com