

Labycare-Tilt Table Manual Adult (SKSS-PP389)



Tilt Table Manual Adult Manufacturer by LABYCARE

Description :

Manually operated Tilt Table with foam padded upholstery, works to standing position manually with easy to turn worm gear arrangement that locks it at any required angle.

Technical :

Range of Tilt : Adjustable from 0 to 90 degree.

Table Top : 195cm long X 61cm wide X 80cm high.

Tilt on Two heavy duty ball bearings.

Fitted with Tilt angle indicator.

Fitted with 35mm thick foam padded upholstery.

Fitted with three, wide straps to hold patient.

Base Frame : Welded frame is 75mm wide for more stability.

Made of rectangular & 30mm square steel tubes.

Legs fitted with 10cm dia. Lockable Castors for Mobility.

Foot Board : Foam Padded removable footboard for easy Entry & Exit.

Worm Gear : Sealed heavy duty system fitted with Four Ball bearings and well Lubricated for life time performance, Works with removable handle.

Finish : Powder coated finish in light cream color.

Weight Capacity : Maximum patient weight 120kg.

We are Exporters and Suppliers of Tilt Table Manual Adult at wholesale competitive prices from India.

Specifications :

Manually operated Tilt Table with foam padded upholstery, works to standing position manually with easy to turn worm gear arrangement that locks it at any required angle.

Technical :

Range of Tilt : Adjustable from 0 to 90 degree.

Table Top : 195cm long X 61cm wide X 80cm high.

Tilt on Two heavy duty ball bearings.

Fitted with Tilt angle indicator.

Fitted with 35mm thick foam padded upholstery.

Fitted with three, wide straps to hold patient.

Base Frame : Welded frame is 75mm wide for more stability.

Made of rectangular & 30mm square steel tubes.

Legs fitted with 10cm dia. Lockable Castors for Mobility.

Foot Board : Foam Padded removable footboard for easy Entry & Exit.

Worm Gear : Sealed heavy duty system fitted with Four Ball bearings and well Lubricated for life time performance, Works with removable handle.

Finish : Powder coated finish in light cream color.

Weight Capacity : Maximum patient weight 120kg.

We are Exporters and Suppliers of Tilt Table Manual Adult at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com