

Labycare-Hi-Low Mat Platform (SKSS-PP391)



Hi-Low Mat Platform Manufacturer by LABYCARE

Description :

HI-LOW Mat Platform is ideal for range of motion, transfer, strength building and other active therapies. Platform with height adjustment and adjustable backrest makes delivering therapy easier for the therapist.

Technical :

Construction : Heavy Steel Frame, with heavy duty '8' bush-bearings for easy operation.

Top : Upholstered, Size 200cm long X 100cm wide. Height Adjustable with Hand Held Switch.

Head end : With Ball Breathing Hole & Plug. Head end elevates for manipulation techniques.

Backrest : Backrest Elevates so that the back can be raised for manipulation techniques & more comfortable sitting.

Height : Table Height is adjustable from 50cm to 95cm. Table provides smooth quite lift with hand held switch.

Lifting System : Fitted with the extremely powerful motor (imported, jack system) that safely accommodates heavy patients (120kg). And ensures rapid and noiseless height adjustments.

Finish : Powder coated frame, beautifully finished.

Power : Unit operates on 220V AC.

We are Exporters and Suppliers of Hi-Low Mat Platform at wholesale competitive prices from India.

Specifications :

HI-LOW Mat Platform is ideal for range of motion, transfer, strength building and other active therapies. Platform with height adjustment and adjustable backrest makes delivering therapy easier for the therapist.

Technical :

Construction : Heavy Steel Frame, with heavy duty '8' bush-bearings for easy operation.

Top : Upholstered, Size 200cm long X 100cm wide. Height Adjustable with Hand Held Switch.

Head end : With Ball Breathing Hole & Plug. Head end elevates for manipulation techniques.

Backrest : Backrest Elevates so that the back can be raised for manipulation techniques & more comfortable sitting.

Height : Table Height is adjustable from 50cm to 95cm. Table provides smooth quite lift with hand held switch.

Lifting System : Fitted with the extremely powerful motor (imported, jack system) that safely accommodates heavy patients (120kg). And ensures rapid and noiseless height adjustments.

Finish : Powder coated frame, beautifully finished.

Power : Unit operates on 220V AC.

We are Exporters and Suppliers of Hi-Low Mat Platform at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com