

Labycare-Stand-In-Frame Child (SKSS-PP396)



Stand-In-Frame Child Manufacturer by LABYCARE

Description :

The durable, functional vertical stander with four points of support provides the therapeutic benefits of standing and places children at peer level. Suitable for tall children.

Features :

Height adjustable chest, pelvic & knee supports.

Adjustable foot restraints.

Detachable rimmed tray for activities.

Chrome plated uprights fitted on non-slip matted foot board.

We are Exporters and Suppliers of Stand-In-Frame Child at wholesale competitive prices from India.

Specifications :

The durable, functional vertical stander with four points of support provides the therapeutic benefits of standing and places children at peer level. Suitable for tall children.

Features :

Height adjustable chest, pelvic & knee supports.

Adjustable foot restraints.

Detachable rimmed tray for activities.

Chrome plated uprights fitted on non-slip matted foot board.

We are Exporters and Suppliers of Stand-In-Frame Child at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com