

Labycare-CP Chair with Activity Tray (SKSS-PP399)



CP Chair with Activity Tray Manufacturer by LABYCARE

Description :

C.P. Rehabilitation Chair with Activity Tray is designed with inclinable foam padded seat & back & fitted with adjustable head support for children who have significant locomotor deficits.

The Chair is most suitable for comfortable positioning of the children with special needs.

Technical :

Frame : Light weight wooden frame with enamel finish.

Seat Tilt : Seat Tilt is adjustable up to 10 degrees with hand levers.

Seat Width : 35cm wide foam padded seat for comfort of child.

Backrest Recline : Adjustable recline upto 25-30 degree for comfortable back support, Useful for the children, who have tendency to bend forward.

Headrest : Foam Padded & Adjustable in Height.

Footrest : Detachable & Adjustable in height. Fitted with feet stabilizing straps & leg support belt.

Wheels : Chair is mounted on Four 5cm dia. Casters with brakes.

Safety : Back is fitted with chest straps to support & hold the child.

Finish : Metal parts have powder coated finish.

We are Exporters and Suppliers of CP Chair with Activity Tray at wholesale competitive prices from India.

Specifications :

C.P. Rehabilitation Chair with Activity Tray is designed with inclinable foam padded seat & back & fitted with adjustable head support for children who have significant locomotor deficits.

The Chair is most suitable for comfortable positioning of the children with special needs.

Technical :

Frame : Light weight wooden frame with enamel finish.

Seat Tilt : Seat Tilt is adjustable up to 10 degrees with hand levers.

Seat Width : 35cm wide foam padded seat for comfort of child.

Backrest Recline : Adjustable recline upto 25-30 degree for comfortable back support, Useful for the children, who have tendency to bend forward.

Headrest : Foam Padded & Adjustable in Height.

Footrest : Detachable & Adjustable in height. Fitted with feet stabilizing straps & leg support belt.

Wheels : Chair is mounted on Four 5cm dia. Casters with brakes.

Safety : Back is fitted with chest straps to support & hold the child.

Finish : Metal parts have powder coated finish.

We are Exporters and Suppliers of CP Chair with Activity Tray at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com