

## Labycare-Finger Goniometer (SKSS-PP405)



Finger Goniometer Manufacturer & Supplier

### **Description :**

A Simple & accurate device for measuring range of motion of small finger joints (I.P. and M.P. joints)  
Constructed of stainless steel.

Comes in a pouch for storage.

Unit measures about 135mm in length.

We are Exporters and Suppliers of Finger Goniometer at wholesale competitive prices from India.

### **Specifications :**

A Simple & accurate device for measuring range of motion of small finger joints (I.P. and M.P. joints)  
Constructed of stainless steel.

Comes in a pouch for storage.

Unit measures about 135mm in length.

We are Exporters and Suppliers of Finger Goniometer at wholesale competitive prices from India.

---

labycare

www.labycare.com | Email: labycare@gmail.com