

Labycare-Gravity Inclinator 360° (SKSS-PP408)



Gravity Inclinator 360° Manufacturer by LABYCARE

Description :

This Goniometer can accurately measure the movement of most of the body joints.

This easy to use Device is placed on the appropriate body part, closet to the axis of movement, with hook and loop strap.

Ideal for measuring typical difficult areas such as the cervical spine (rotation, flexion, extension & lateral flexion) and lumbar spine.

We are Exporters and Suppliers of Gravity Inclinator 360° at wholesale competitive prices from India.

Specifications :

This Goniometer can accurately measure the movement of most of the body joints.

This easy to use Device is placed on the appropriate body part, closet to the axis of movement, with hook and loop strap.

Ideal for measuring typical difficult areas such as the cervical spine (rotation, flexion, extension & lateral flexion) and lumbar spine.

We are Exporters and Suppliers of Gravity Inclinator 360° at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com