

Labycare-Upper Extremity Exercise Station (SKSS-PP410)



Upper Extremity Exercise Station Manufactruer by LABYCARE

Description :

The Sturdy, mobile exercise station contains everything you need for optimum patient therapy and provides storage space for the following fourteen most recognized activities appropriate to a variety of patients to improve : ROM, Flexibility, Strengthening and Endurance of fingers, grip, wrist, elbow & shoulder.

Contents :

- * **Dumbbells Set** : One pair each of 1 kg, 1.5kg, 2kg, 2.5kg & 3kg. (Total Weight 20kg)
- * **Cuff Weights** : One pair each of 0.5kg, 1kg, 1.5kg, 2kg & 3kg (Weight 16kg)
- * **Resistive Bands** : Five Resistance in 5.5 meter length.
- * **Therapy Putty** : Five Packs of 100gm each in Five different resistance.
- * **Grip Exerciser** : One, with six springs.
- * **Medicine Ball Set** : Five balls, 1,2,3,4 & 5kg.
- * **Hand Gym Kit Board** : One, for fingers & grip exercises.
- * **Wrist Roll with Supination & Pronation** : One
- * **Hand Skate Roller** : One
- * **Shoulder Pulley Kit** : One
- * **Spring with Marked Tension** : Set of Five Springs.

* **Shoulder Abduction Ladder** : One

* **Wire Handle** : Two

* **Holding Mitt** : One

Construction : Hardwood & Ply board construction with three storage shelves & mounted on four 4" dia. Lockable casters. **Approx. Size** : 28" x 28"x 68" high.

We are Exporters and Suppliers of Upper Extremity Exercise Station at wholesale competitive prices from India.

Specifications :

The Sturdy, mobile exercise station contains everything you need for optimum patient therapy and provides storage space for the following fourteen most recognized activities appropriate to a variety of patients to improve : ROM, Flexibility, Strengthening and Endurance of fingers, grip, wrist, elbow & shoulder.

Contents :

* **Dumbbells Set** : One pair each of 1 kg, 1.5kg, 2kg, 2.5kg & 3kg. (Total Weight 20kg)

* **Cuff Weights** : One pair each of 0.5kg, 1kg, 1.5kg, 2kg & 3kg (Weight 16kg)

* **Resistive Bands** : Five Resistance in 5.5 meter length.

* **Therapy Putty** : Five Packs of 100gm each in Five different resistance.

* **Grip Exerciser** : One, with six springs.

* **Medicine Ball Set** : Five balls, 1,2,3,4 & 5kg.

* **Hand Gym Kit Board** : One, for fingers & grip exercises.

* **Wrist Roll with Supination & Pronation** : One

* **Hand Skate Roller** : One

* **Shoulder Pulley Kit** : One

* **Spring with Marked Tension** : Set of Five Springs.

* **Shoulder Abduction Ladder** : One

* **Wire Handle** : Two

* **Holding Mitt** : One

Construction : Hardwood & Ply board construction with three storage shelves & mounted on four 4" dia. Lockable casters. **Approx. Size** : 28" x 28"x 68" high.

We are Exporters and Suppliers of Upper Extremity Exercise Station at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com