

Labycare-Arm Sling Strap With Shoulder Cushion (SKSS-RP434)



Arm Sling Strap With Shoulder Cushion Manufacturer by LABYCARE

Description :

Arm Sling Strap helps in management of strain, sprain and fracture.

Arm Sling Strap supports and suspends the forearm with elbow in 90° flexion.

Adjustable straps according to need and neck pad for even pressure distribution.

Easy to wear.

Arm Sling Strap is common for left & right hand.

We are Exporters and Suppliers of Arm Sling Strap With Shoulder Cushion at wholesale competitive prices from India.

Specifications :

Arm Sling Strap helps in management of strain, sprain and fracture.

Arm Sling Strap supports and suspends the forearm with elbow in 90° flexion.

Adjustable straps according to need and neck pad for even pressure distribution.

Easy to wear.

Arm Sling Strap is common for left & right hand.

We are Exporters and Suppliers of Arm Sling Strap With Shoulder Cushion at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com