

Labycare-Contoured L S Support (SKSS-RP444)



Contoured L S Support Manufacturer by LABYCARE

Description :

Contoured L S Support Double elastic strap design with four contoured malleable aluminium splints.

Contoured L S Support has four aluminium malleable metal splints which are removable & can be given the desired contour to support the spine.

Contoured L S Support helps in low back pain, post surgery in lumbosacral region and other spinal conditions.

Contoured L S Support ensures better immobilization with good posture.

Contoured L S Support boarder elastic gives near rigid immobilization.

We are Exporters and Suppliers of Contoured L S Support at wholesale competitive prices from India.

Specifications :

Contoured L S Support Double elastic strap design with four contoured malleable aluminium splints.

Contoured L S Support has four aluminium malleable metal splints which are removable & can be given the desired contour to support the spine.

Contoured L S Support helps in low back pain, post surgery in lumbosacral region and other spinal conditions.

Contoured L S Support ensures better immobilization with good posture.

Contoured L S Support boarder elastic gives near rigid immobilization.

We are Exporters and Suppliers of Contoured L S Support at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com