

Labycare-Dorsolumbar Spinal Brace (SKSS-RP447)



Dorsolumbar Spinal Brace Manufacturer by LABYCARE

Description :

Dorsolumbar Spinal Brace is recommended in vertebral body fractures or collapse, osteoporosis, spondylolisthesis, post-up spinal fusion, laminectomy and disectomy. Supports and immobilizes spine in neutral position.

Dorsolumbar Spinal Brace is made of posterior contoured and malleable metal frame with pelvic band, thoracic band and two paraspinous uprights.

Dorsolumbar Spinal Brace has abdominal attachment and shoulder straps to hold the brace snugly. Shoulder straps provides proper fitting and supports immobilization.

Dorsolumbar Spinal Brace is ideal for weakened dorsolumbar spine due to TB & suitable for post operative rehabilitation of back ache patient.

We are Exporters and Suppliers of Dorsolumbar Spinal Brace at wholesale competitive prices from India.

Specifications :

Dorsolumbar Spinal Brace is recommended in vertebral body fractures or collapse, osteoporosis, spondylolisthesis, post-up spinal fusion, laminectomy and disectomy. Supports and immobilizes spine in neutral position.

Dorsolumbar Spinal Brace is made of posterior contoured and malleable metal frame with pelvic band, thoracic band and two paraspinous uprights.

Dorsolumbar Spinal Brace has abdominal attachment and shoulder straps to hold the brace snugly. Shoulder straps provides proper fitting and supports immobilization.

Dorsolumbar Spinal Brace is ideal for weakened dorsolumbar spine due to TB & suitable for post

operative rehabilitation of back ache patient.

We are Exporters and Suppliers of Dorsolumbar Spinal Brace at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com