

Labycare-philadelphia collar (SKSS-RP450)



philadelphia collar Manufacturer by LABYCARE

Description :

Philadelphia Collar is also known as Cervical Orthosis

Philadelphia Collar is used during fractures of cervical vertebrae, post trauma, disc herniation, spondylolisthesis and other similar conditions.

Philadelphia Collar supports and immobilizes the neck region with comfort 2 piece design with anterior and posterior struts of semi rigid plastic.

Philadelphia Collar has molded chin and occipital support. Offers large trachea opening.

We are Exporters and Suppliers of Philadelphia Collar or Cervical Orthosis at wholesale competitive prices from India.

Specifications :

Philadelphia Collar is also known as Cervical Orthosis

Philadelphia Collar is used during fractures of cervical vertebrae, post trauma, disc herniation, spondylolisthesis and other similar conditions.

Philadelphia Collar supports and immobilizes the neck region with comfort 2 piece design with anterior and posterior struts of semi rigid plastic.

Philadelphia Collar has molded chin and occipital support. Offers large trachea opening.

We are Exporters and Suppliers of Philadelphia Collar or Cervical Orthosis at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com