

Labycare-Kymograph Manual (SKSSKM912B)



Kymograph Manual (Analog)

Description :

Kymograph Manual (Analog)

Specifications :

Kymograph (Analog)	Manual	ELECTRICAL	E-8	SHERRINGTON
		KYMAGRAPH; oil lapped brass gear system; speeds from 0.12mm/ sec. to 640 mm/sec. in 8 steps. Jerk free running, 8 gear changes through a dial type knob; stainless steel spindle with screw lift. 15x15 cm. cylinder with pilot lamp, plug & cord		
TECHNICAL SPECIFICATION:	• Wide speed variation - from 1 rev. in 25 min to 1 rev in 11 hours. • Permissible paper length 2.5 meters, adjustable along precision machined 120 cm. Bed. • Automatic paper tensioning devices & flat writing surface attachment. • Neutralisation clutch for instantaneous start or stop. • Centralised electrical points for Respiration Pump. Time Clock & Operation Table. • Mounted on Heavy cast metal platform with castors for mobility. • With Manometer & two way Time Marker.			

labycare

www.labycare.com | Email: labycare@gmail.com