

# Labycare-Patient Care Mattress (SKSS-809J)



Patient Care Mattress Manufacturer by LABYCARE

## **Description :**

The ICU bed mattress is designed for optimal comfort and hygiene in critical care settings. With dimensions of 1970 mm (L) x 860 mm (W) x 100 mm (T) and a core made of high-density PU foam (40 density), it ensures durability and support for patients. The mattress features a waterproof, fire-retardant, and anti-microbial PU/PVC-coated fabric cover that is easy to clean, enhancing infection control. Its four-section design perfectly accommodates fowler beds, allowing smooth adjustments for backrest and knee rest positions. Supporting a load capacity of up to 150 kg, this ISO 13485-certified mattress meets stringent quality and safety standards, making it ideal for ICU environments.

## **Specifications :**

**Dimensions:** 1970 mm (L) x 860 mm (W) x 100 mm (T).

**Material:** High-density PU foam with 40 density for durability and patient support.

**Cover Fabric:** PU/PVC-coated, waterproof, fire-retardant, and anti-microbial for hygiene and easy cleaning.

**Design:** Four-section construction, compatible with fowler beds for seamless backrest and knee rest adjustments.

**Load Capacity:** Supports up to 150 kg.

**Compliance:** ISO 13485-certified, ensuring high safety and quality standards.

---

labycare

www.labycare.com | Email: labycare@gmail.com