

Labycare-Hyper Extension Brace (SKSS-RP460)



Hyper Extension Brace Manufacturer by LABYCARE

Description :

Hyper Extension Brace is light weight adjustable, comfortable & easy to apply with Velcro strap, sliding arrangement.

Hyper Extension Brace consist of a rectangular metal frame that goes over the front of the body & puts pressure on the chest and pubic bone & keeps the spine in an extended position.

Hyper Extension Brace with 3 pads help to stabilize the spine and prevent forward movement. One pad is located along the abdominals, another is higher up on the chest & the third pad is on the back & covers the affected area & restricts forward bending.

Hyper Extension Brace is recommended in undisplaced stable fractures and (compression) fractures of the Lower thoracic & upper lumbar vertebrae, any complication of thoracolumbar junction which needs moderate immobilization disc disease, to prevent deformity in Ankylosing spondylitis, Post Operative spine immobilization etc.

We are Exporters and Suppliers of Hyper Extension Brace at wholesale competitive prices from India.

Specifications :

Hyper Extension Brace is light weight adjustable, comfortable & easy to apply with Velcro strap, sliding arrangement.

Hyper Extension Brace consist of a rectangular metal frame that goes over the front of the body & puts pressure on the chest and pubic bone & keeps the spine in an extended position.

Hyper Extension Brace with 3 pads help to stabilize the spine and prevent forward movement. One pad is located along the abdominals, another is higher up on the chest & the third pad is on the back & covers the affected area & restricts forward bending.

Hyper Extension Brace is recommended in undisplaced stable fractures and (compression) fractures of the Lower thoracic & upper lumbar vertebrae, any complication of thoracolumbar junction which needs moderate immobilization disc disease, to prevent deformity in Ankylosing spondylitis, Post Operative spine immobilization etc.

We are Exporters and Suppliers of Hyper Extension Brace at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com