# Labycare-Vaporizer Halothane Standard (SKSS-AP133)



Vaporizer Halothane Standard Manufacturer By LABYCARE

**Description :**

Vaporizer Halothane can only fill with Halothane

Range of concentration Halothane (0.5-5) vol %

Concentration error: ±0.25vol% or 20% (the larger value)

Range of flow compensation: 0.5 L/min~10L/min

Range of temperature compensation: 15 degree celcius~35 degree celcius

The greatest angle of slope is less than 45 degrees

Injection volume: the lowest scale mark is about 30mL and the highest scale mark is 120mL

Dimensions : Height 225mm, Width 105 mm, Depth 210 mm.

We are Manufacturer, Exporters and Suppliers of Vaporizer Halothane Standard at wholesale competitive prices from India.

**Specifications :**

Vaporizer Halothane can only fill with Halothane

Range of concentration Halothane (0.5-5) vol %

Concentration error: ±0.25vol% or 20% (the larger value)

Range of flow compensation: 0.5 L/min~10L/min

Range of temperature compensation: 15 degree celcius~35 degree celcius

The greatest angle of slope is less than 45 degrees

Injection volume: the lowest scale mark is about 30mL and the highest scale mark is 120mL

Dimensions : Height 225mm, Width 105 mm, Depth 210 mm.

We are Manufacturer, Exporters and Suppliers of Vaporizer Halothane Standard at wholesale competitive prices from India.

---