Labycare-scrotal support (SKSS-RP420)



scrotal support Manufacturer By LABYCARE

Description :

Scrotal Support is recommended following surgery and in treatment of many scrotal problems.

Scrotal Support helps to provide support to the Scrotum cremaster muscle which adds comfort and convenience for the patient.

Scrotal Support is light adjustable waistband and soft knit stretchable pouch of polyester. It is made of soft and hypoallergenic material.

Specifications :

Scrotal Support is recommended following surgery and in treatment of many scrotal problems.

Scrotal Support helps to provide support to the Scrotum cremaster muscle which adds comfort and convenience for the patient.

Scrotal Support is light adjustable waistband and soft knit stretchable pouch of polyester. It is made of soft and hypoallergenic material.